

**D-GLUCOSE
(LIQUID READY™)
VALIDATION REPORT**

**SKU: 700007404
K-GLULQR**

11/24



INTRODUCTION:

The D-Glucose Assay Kit (Liquid Ready™) is a robust, quick and easy method for the measurement of D-Glucose in various matrices and is fully automatable for high throughput analysis of samples. Data presented in this report validates that this method is fit for the purpose intended.

RECOMMENDATIONS FOR ANALYSIS:

Please reach out to your local sales representative or to the technical team should you require any assistance, particularly in relation to assay troubleshooting, data analysis, additional matrix testing and application support in relation to automated analysers.

- This test should only be carried out by trained laboratory employees. The product instructions must be followed to ensure an accurate and robust result.
- Store the kit at 2 - 8 °C.
- Ensure all kit components come to room temperature 20 - 25 °C before use.
- Use the contents of bottles 1, 2 & 3 as supplied.
- Use a repetitive pipettor to add reagents to the cuvette for increased accuracy.
- The reagent blank value must be determined once for each run and subtracted from each sample result.
- Users should perform matrix validation experiments prior to routine use. This process will highlight any problematic matrices encountered.
- The Carrez Clarification kit (K-CARREZ; 700004270) is recommended for fat removal and deproteinization if necessary.
- Use separate tips for each sample extract and control solutions to avoid cross-contamination. Additionally, pre-flush the tip before pipetting.
- When testing solid samples, ensure a representative portion is taken and homogenized before weighing.

EQUIPMENT (RECOMMENDED):

Matrix Extraction

1. Syringe filters (AGILENT™ Nylon or equivalent, 0.2 micron); alternatively, eppendorf tubes and centrifuge.
2. Sartorius™ grade 292 filter paper or equivalent.
3. Blender (e.g. Nutribullet™).
4. Analytical balance and spatulas.
5. Wide-mouth bottles (e.g. Duran® bottles, 100 mL and 1 L).
6. Glass graduated cylinders (25 mL, 50 mL) and volumetric flasks (100 mL).
7. pH-meter or equivalent.
8. Refrigerator or ice bath (or equivalent).
9. Stir plate.

SUMMARY OF PERFORMANCE DATA:


All testing was performed using the standard manual assay described in the product instruction document. Results are summarised in the table below:

Recommended working range	0.010 g/L – 1.500 g/L
Limit of Detection (LOD)	0.002 g/L
Limit of Quantification (LOQ)	0.006 g/L
Limit of Precision	0.010 g/L
Specificity	D-Glucose
Bias %	-1.2% to 1.5%
Acceptable Recovery of Standards	95% - 105%
Stability Studies	24 months shelf life from date of manufacture, see product label for expiry
	Kit performance maintained after 3 freeze-thaw cycles
Robustness	< 5 mins to reach completion at 20°C, 25°C and 37°C
Repeatability	CV < 5% for pure D-Glucose samples
	CV < 10% for a range of matrices tested
Selectivity/Cross-reactivity	Cross reactivity observed with D-Fructose above 100 mg/mL in assay
Matrix Interference	Recovery between 95% - 105% in red wine, white wine, apple juice, orange juice, sour jellies, fruit liquid yogurt, grapes and breakfast cereals

LINEARITY AND WORKING RANGE:

The recommended linear measurement range is from 0.01 g/L to 1.5 g/L of D-Glucose. Samples with concentrations up to 1.75 g/L can be analyzed without dilution, however assay repeatability/robustness may be affected at concentrations above 1.5 g/L. Samples containing higher levels of analyte can be diluted with distilled water before measurement. Values below 0.01 g/L will not fall within the acceptance criteria.

D-Glucose (g/L)	$\Delta_{\text{absorbance}}$ (340 nm)	Recovery (%)
0	0	0
0.006	0.008	110.2
0.010	0.010	100.37
0.050	0.051	101.12
0.100	0.101	101.05
0.400	0.393	98.162
0.800	0.816	102.00
1.000	1.008	100.81
1.500	1.500	99.973
1.750	1.706	97.486
2.000	1.851	92.527



LIMIT OF DETECTION, QUANTIFICATION AND PRECISION:

The LOD is the lowest concentration of the analyte that can be detected by the method. This was determined by testing 20 replicates of the blank (i.e. adding 100 μL of water instead of sample). The $\Delta A_{\text{Limit of Detection}}$ is calculated as $3.3 \times s'0$; where $s'0$ is the standard deviation of a number of samples $\Delta_{\text{Absorbance}}$ reading.

The LOQ is the lowest level at which the kit's performance is acceptably repeatable. This was determined by testing 20 replicates of the blank (i.e. adding 100 μL of water instead of sample). The $\Delta A_{\text{Limit of Quantification}}$ is calculated as $kQ \times s'0$; where $s'0$ is the standard deviation of a number of samples $\Delta_{\text{Absorbance}}$ reading. The IUPAC default value for kQ is 10.

The Limit of Precision (g/L) is the analyte concentration at which it was experimentally determined that acceptable recoveries ($\pm 5\%$) are routinely achieved. It was assessed by measuring 12 replicates of the lowest analyte concentrations, starting at LOQ level.

$\Delta A_{\text{Limit of Detection}}$	Limit of Detection (g/L)	$\Delta A_{\text{Limit of Quantification}}$	Limit of Quantification (g/L)	Precision Limit (g/L)
0.0025	0.002	0.0075	0.006	0.010

NOTE: The above detection limits were calculated based on assay concentration (i.e. samples post-extraction). The dilution used in pre-treatment must be accounted for when establishing the detection limits for specific samples.

TRUENESS AND BIAS:

The Trueness of the D-Glucose Assay kit (Liquid Ready™) was tested using validated aqueous D-Glucose standards. The *bias* is the comparison of the mean of the results (X) achieved using the standard manual protocol with a suitable reference material.

- Relative Bias is calculated in percent as: $b(\%) = (X - X_{ref}) / X_{ref} \times 100$

Reference material (g/L)	Replicates, n	% CV	% Recovery	% Bias
0.010	12	3.136	101.53	1.527
0.750	12	0.669	100.57	0.570
1.500	12	0.385	98.787	-1.213

Recovery of sample solution was within the acceptance criteria of $100 \pm 5\%$ and %CV below 5%. The *bias* for this method ranges from -1.21 to 1.52%.

INTERFERENCE AND SELECTIVITY:

Selectivity

The assessment of the selectivity of this method towards D-Glucose in the presence of interfering compounds was tested by spiking a set concentration of interfering agent with a known concentration of validated aqueous D-Glucose standard and assessing the recovery. Concentrations of interfering agent were selected based on what were considered to be levels likely to be found in relevant samples.

NOTE: Interfering substances in the sample being analysed can be identified by including an internal standard. Quantitative recovery of this standard would be expected. Losses in sample handling and extraction are identified by performing recovery experiments, i.e. by adding D-Glucose to the sample in the initial extraction steps.

Compound tested	Concentration of interferent in assay (g/L)	Recovery of D-Glucose (%)
D-Fructose	150	86.15
	100	98.75
D-Galactose	10	99.60
D-Mannose	5	99.7
Sucrose	200	100.2
Lactose	20	101.7
Glycerol	10	100.1
Sorbitol	200	100.1
L-Malic acid	10	100.0
L-Tartaric acid	10	99.6
Citric acid	5	100.0
Maleic Acid	5	100.2
Succinic Acid	5	100.5
Ascorbic acid	50	99.5
Galacturonic acid	10	100.4
D/L-Lactic acid	20	100.4
Acetic acid	50	98.9
Ethanol	120	99.7
Sodium Chloride	20	99.5
Calcium Chloride	100	99.3

Most recoveries were found to be within the acceptance criteria of $100 \pm 5\%$. The analysis shows that there is some

cross reactivity with Fructose above 100 g/L in assay. This means that samples containing high levels of fructose will need to be diluted prior assay. No other cross-reactivity was found in any of the tested substances.

Matrix Interference

The assessment of the interference of this method when using specific complex matrices was tested by spiking extracted samples with a known concentration of validated aqueous ammonia standard and assessing the recovery.

NOTE: All matrices were prepared according to the sample preparation methods described in the product instruction document which can be found on the product webpage.

Matrix tested	Replicates, n	D-Glucose Spike Recovery (%)
Red Wine	2	102.3
White Wine	2	100.3
Apple Juice	2	102.3
Orange Juice	2	100.8
Breakfast Cereals	2	98.97
Sour Jellies	2	100.03
Fruit Yogurt	2	101.61
Grapes	2	100.47

Eight matrices were tested; all relevant samples fall within 95-105% recovery. Furthermore, the average %CV for all spike recoveries was 2.28% which demonstrates the precision and specificity of this method.

ROBUSTNESS AND STABILITY

Storage Temperature

To assess the storage stability of the test kit components, all kit components were stored at 4°C. Real time monthly performance and enzyme activity testing was conducted, and slope predictions were used to estimate the product shelf life.

Storage Temperature	Reagent Tested	Stability Data
4°C	Reagent 1	24 months shelf life from date of manufacture, see product label for expiry
	Reagent 2	

The storage robustness of this kit was also tested by performing three freeze-thaw cycles. All kit components were frozen at -20°C overnight and allowed to thaw before testing. This freeze-thaw cycle was repeated 3 times using the same components.

D-Glucose Recovery (g/L)				
Expected D-Glucose (g/L)	T0	1st Cycle	2nd Cycle	3rd Cycle
0.010	0.010	0.010	0.010	0.010
0.750	0.746	0.739	0.738	0.732
1.500	1.477	1.456	1.484	1.464

Recovery of sample solutions was within the acceptance criteria of $100 \pm 5\%$. ANOVA evaluation of the freeze-thaw cycles shows that there is no significant difference in assay performance over three consecutive freezing cycles of K-GLULQR reagents.

Assay Temperature

Enzymatic test kits are known to be subject to variation when tested under different environmental conditions; temperature being a key parameter that influences reaction rate and analyte recovery. The D-Glucose Assay kit (Liquid Ready™) was tested using different reaction temperatures (20°C, 25°C and 37°C); recoveries and reaction times were analysed.

Expected D-Glucose (g/L)	D-Glucose Recovery (g/L)		
	20°C	25°C	37°C
0.010	0.010	0.010	0.010
0.750	0.739	0.746	0.740
1.500	1.482	1.477	1.478
Reaction Time	4 min	4 min	2 min
<i>Recommended time</i>	<i>5 minute reaction time is recommended</i>		

ANOVA evaluation of the effect of temperature in assay shows that there is no significant difference in recovery at the temperatures tested. All recoveries were found to be within acceptance criteria of $100 \pm 5\%$ with a faster reaction time observed at 37°C. The recommended reaction time for K-GLULQR is 5 minutes.

PRECISION AND REPEATABILITY

Precision

Precision is a measure of the variability in results, on different days and by different analysts over a period of time and using different lots of the test kit. The precision of the D-Glucose Assay kit (Liquid Ready™) was investigated using validated aqueous D-Glucose samples.

D-Glucose (g/L)	Replicates, n	Mean (g/L)	Standard Deviation	%CV
0.010	24	0.010	0.0003	2.693
0.750	24	0.749	0.0032	0.433
1.500	24	1.496	0.006	0.398

All recoveries were shown to be within acceptance criteria of $100 \pm 5\%$ and %CV below 5%. This demonstrates that the method is precise and repeatable when using validated aqueous D-Glucose standards.

NOTE: Users should perform matrix validation work prior to routine use. This process will highlight any problematic matrices encountered. If you have questions about these or other matrices, please contact your local sales representative for support.

Matrix Repeatability

Matrix repeatability was tested by one analyst over three days using a range of selected matrices. Sample preparation was performed on a daily basis to assess precision of analyte extraction as well as precision of the test method. All extractions were performed as per the sample preparation method stated in the product instruction.

Matrix	Replicates, n	Mean result (g/L or g/100g)	Standard Deviation	%CV
Red wine	6	0.354	0.031	8.750
White wine	6	0.820	0.015	1.883
Apple juice	6	20.46	1.355	6.619
Breakfast Cereal	6	0.614	0.007	1.081
Orange Juice	6	18.73	0.533	2.847
Sour Jellies	6	11.60	0.357	3.074
Fruit Yogurt	6	1.537	0.075	4.863
Grapes	6	7.654	0.364	4.756


All sample analysis resulted in %CV values of <10%. This demonstrates that the method is precise and repeatable when testing complex matrices.

METHOD AUTOMATION

The D-Glucose Liquid Ready™ kit was designed for biochemistry analyzers. It is easily adapted to any instrument. D-Glucose content can be determined in a single test using a cubic spline fit.

NOTE: For each batch of samples that is applied to the determination of D-Glucose a calibration curve must be performed concurrently using the same batch of reagents.

D-Glucose (g/L)	Replicates, n	Mean result (g/L)	Standard Deviation	%CV
0.015	4	0.151	0.0007	0.410
0.750	4	0.748	0.0333	0.388
1.500	4	1.514	0.0035	0.207



SERVICES AND TECHNICAL SUPPORT

Please reach out to your local sales representative should you require any assistance, particularly in relation to:

- Troubleshooting
- Data analysis
- Additional matrix testing
- Application support in relation to automated analyzers.

Supporting documents can be found in the product page:

- Product Instructions
- D-Glucose Liquid Ready™ Quick Reference Guide
- Mega-Calc™ D-Glucose Liquid Ready™
- Safety data sheets (SDS)
- Certificate of analysis (CoA)



Contact us for more information: neogen.com/contact

Without guarantee

The information contained in this report is, to the best of our knowledge, true and accurate, but since the conditions of use are beyond our control, no warranty is given or is implied in respect of any recommendation or suggestions which may be made or that any use will not infringe any patents. It is the user's responsibility to perform in-house matrix validation work prior to routine use.

© 2024, Neogen Corporation. All rights reserved.

Neogen, Megazyme, Liquid Ready, and Mega-Calc are trademarks of Neogen Corporation. All other trademarks are property of their respective owners.